



UROTECH

Leading Urotechnology Made in Germany

Stop the stones

URETERAL OCCLUSION CATHETER

- Hydrophilic tip for easy positioning
- Additional radiomarker in the centre of the balloon for safe localisation
- Smooth junction of catheter tube and balloon
- Latex free balloon
- Active deflation of balloon with syringe possible
- Easy adjustment of connector

Size Catheter (FR/CH)	Size Tip (FR/CH)	Catheter Length (cm)	Balloon Volume (ml)	Mandril Length (cm)	Mandril Diameter (inch)	Order Code
5.0	4.0	75	1	84	0.035"	OK-05-75