



The solution for confidence, comfort and efficiency in continence care

There are numerous reasons why people need a transurethral or suprapubic urine drainage system for a short or long period. Regardless of whether in clinical use, in geriatric and nursing care or for self-catheterisation in a domestic setting: Dealing with catheters and urine bags is a sensitive issue for everyone involved. On the one hand, there is the unavoidable invasion of the private sphere by professional personnel. On the other hand, there is the fear of stigmatisation due to the loss of bladder function and, of course, the risk of infection, urinary retention, urethral stricture or the leakage of urine and odour. Discretion and confidence through high-quality, comfortable products that are easy to use in every individual case – that's where URELIA[®] comes in. The urine bag brand from the urology specialist UROVISION offers a complete range of products of excellent quality that has been thought through down to the last detail. The cost? Unbelievable.

URELIA[®] – we focus on bladder care so that you can forget about it

- Perfectly tailored to the individual requirements of patients
- High-quality products for holistic care
- Increased wearing comfort due to skin-friendly coating on hoses
- Latex-free
- CE marking

There are many good reasons for taking a closer look at the URELIA[®] range.
See for yourself.

URELIA® by UROVISION:

Urine drainage can be this easy and secure

URELIA® offers the complete range of products for the management of incontinence for every patient's needs and situation. Reliable, practical and all from one source.

Closed urine drainage system for long-term care



Properties:

- └ 2,000 ml volume
- └ Hose length 120 cm, kink resistant with sliding clamp
- └ Needle-free urine sampling site
- └ Drip chamber with integrated non-return valve
- └ Bacteria-proof discharge valve
- └ Simple one-hand drainage
- └ Double hook attachment
- └ Bed sheet clamp
- └ Sterile

Order no.	Volume	Hose length	Unit
UVD-23	2,000 ml	120 cm	20

Urine drainage system for short-term use



Properties:

- └ 2,000 ml volume
- └ Hose length 110 cm, kink resistant with sliding clamp
- └ Needle-free urine sampling site
- └ Return barrier
- └ Bacteria-proof discharge valve
- └ Simple one-hand drainage
- └ Conical step connector for optimum hold
- └ Sterile

Order no.	Volume	Hose length	Unit
UVD-43	2,000 ml	110 cm	50

3-chamber leg bag for high wearing comfort with inconspicuous fit



Properties:

- └ Stable urine bags with 500 ml or 750 ml volume
- └ Adapts inconspicuously to the shape of the leg thanks to the three chambers
- └ Special, skin-friendly, soft back of the bag
- └ Coated hose to avoid irritation
- └ Hose length 50 cm, can be shortened with adapter
- └ Simple one-hand drainage
- └ Return valve
- └ Can be connected to bed bag
- └ Washable special fastening straps for secure attachment to the leg (available separately)
- └ Sterile

Order no.	Volume	Hose length	Unit
UVD-3K50	500 ml	50 cm (to shorten)	50
UVD-3K75	750 ml	50 cm (to shorten)	50

Wide range of high-quality, excellent leg bags with dog house



Properties:

- └ Stable leg bags with 500 ml or 750 ml volume
- └ Hose length 10 cm or 50 cm (can be shortened with exchangeable universal adapter)
- └ Coated hose to avoid irritation
- └ Special, skin-friendly, soft back of the bag
- └ Dog house for easy storage of the drain tap
- └ Washable special fastening straps for secure attachment on the leg (available separately)
- └ Can be connected to bed bag
- └ Simple one-hand drainage
- └ Return valve
- └ Sterile

Order no.	Volume	Hose length	Unit
UVD-5010	500 ml	10 cm	50
UVD-5050	500 ml	50 cm (to shorten)	50
UVD-7510	750 ml	10 cm	50
UVD-7550	750 ml	50 cm (to shorten)	50

The secure all-round care for day and night



Properties:

- └ Sturdy day/night bags with 1,000 ml volume
- └ Special, skin-friendly, soft back of the bag
- └ Coated hose to avoid irritation
- └ Different hose lengths available
- └ Needle-free urine sampling site
- └ Drip chamber with integrated non-return valve
- └ Bacteria-proof discharge valve
- └ Dog house for easy storage of the drain tap
- └ Simple one-hand drainage
- └ Sterile

Order no.	Volume	Hose length	Unit
UVD-1030	1,000 ml	30 cm	50
UVD-1060	1,000 ml	60 cm	50
UVD-1090	1,000 ml	90 cm	50

Fastening straps for leg bags (available separately)

Properties:

- └ Set of two straps
- └ 5.5 cm wide strap and 2.5 cm wide strap
- └ With silicone coating inside for optimum grip
- └ Washable at 30° Celsius

Order no.	Unit
UVD-552500	100





Stay on the move

Worldwide around 200 million people suffer from urinary incontinence for various reasons. If an absorbent treatment is not sufficient, transurethral catheters and – especially in the case of neural micturition disorder – suprapubic catheters connected to a urine collection bag are used. A urine drainage system must support the mobility of independent people and remain unnoticeable. This is achieved by using high-quality materials that offer a high level of comfort and a feeling of security. Essential aspects for those affected: Discretion and confidence while preserving quality of life and mobility.

With superbly finished URELIA® products you are in good hands, thanks to high-quality materials that are kind to your skin and have anti-adherent surfaces. URELIA® is the urine bag brand of UROVISION. UROVISION offers a wide range of transurethral and suprapubic catheters that have perfect compatibility and outstanding material quality.

Enabling patients to live a dynamic and agile life – even if your bladder lets you down: This is where we come in.

Please also contact us for suitable catheters for the URELIA® urine bags.



UROVISION – UROTECH

Mastering medical challenges by innovation.

URELIA® by UROVISION

Urovision Gesellschaft für medizinischen Technologie-Transfer mbH

Medi-Globe-Straße 1-5

D-83101 Rohrdorf OT Achenmühle

Phone: +49 (0) 80 32 973 - 400

Fax: +49 (0) 80 32 973 - 411

E-Mail: urovision@urovision.de

www.urovision.de

Art.-No.: LL-0239 UV PROMO Urelia_EN
Vers. 1.0 11/2017