



UROTECH

Leading Urotechnology Made in Germany

No barrier for SEMISLIP

SEMISLIP GUIDE WIRE suitable for all urological procedures

- 9 cm, very flexible and hydrophilic tip provides none-traumatic insertion and safe bypass of obstructions
- PTFE-coated shaft supports easy and fast placement and reduces the risk of dislocation during urological procedures

Order Codes	Diameter (inch)	Length (cm)	Tip	Tip Stiffness	Shaft Stiffness
UAS-35150-N	0.035"	150	straight	standard	soft
UAS-35150-S	0.035"	150	straight	standard	standard
UAS-35150-SS	0.035"	150	straight	standard	stiff